

# SonoScape

## E3

A brand new portable ultrasound system E3, brings you a distinct experience with SonoScape's traditional imaging technologies.

E3's accurate B mode and sensitive color signal give crisp, detailed images to improve your scanning experience while increasing your diagnostic confidence.



# Excellent Image Quality as Always

## C-Field Beam

The continuously dynamic focus provides more energy which contributes to higher contrast resolution, signal-noise ratio and uniformity in the image.

## $\mu$ -Scan

The latest generation of  $\mu$ -Scan imaging greatly enhances the image by reducing noise, improving signal strength and improving visualization.

## Dynamic Color

With Dynamic Color the sonographers can easily see in detail very small veins and slower velocities for detailed blood flow information of the patient.

## Tissue Specific Imaging

The system detects different tissue automatically by matching different acoustic ranges, from which the user can then acquire images with more uniformity and higher spatial resolution.

# Wide Range of Probes and Applications

