




The companion animal package provides users with complete solutions for a variety of applications and animals. From the high definition screen and control panel design to the boot up speed and preset examination options, the X3V accessories the examination process.

Efficient Experience

Since X3V adopts SonoScape's latest imaging technologies to produce excellent images with high resolution and deep penetration, presenting a large amount of information onto the ultrasound screen. Improved image quality will allow for faster and more accurate exams of any species.

- Multi-beam processing technology
- C-field beam technology
- Spatial compound imaging
- μ -Scan image Processing technology
- Pure inversion harmonic imaging



15.6-inch high definition monitor with large display region



Noiseless system to prevent startling the animals




Boot up quickly with predetermined exam settings loaded and ready



Auto brightness adjustment screen and backlit control panel



Keyboard cover to avoid fluids and dust



SonoScape

SonoScape Medical Corp.
2F, 12th Building, Shenzhen Software Park Phase II, Keji Middle
2nd Road, Shenzhen, China
Tel: +86-755-26722890 Fax: +86-755-26722850
Email: market@sonoscape.net www.sonoscape.com

SonoScape

X3V

