

Balanced Design Precise Image
When Every Pawprint Counts

ProPet X10/E10

VET



Powered by C-Field⁺ Platform

Delivering unmatched image clarity and ultra-fast signal processing, it ensures precise, real-time diagnostics you can rely on.



Dedicated Veterinary Software for Animal Care

A fully integrated software tailored for veterinary imaging, clinical workflow, and patient data management, featuring 200+ veterinary bodymarks, 120+ annotations, and a comprehensive library of 20+ species.

Flexible, Intuitive & Customizable Veterinary Exam Mode



Measure



Report



Species-Specific
BSA Formulas

Integrated Measurement, Report & BSA Tools for Veterinary Applications



SonoHelp

Designed specifically for veterinarians, this embedded software provides reference imaging planes and guidance on proper scanning techniques, with support for custom additions of images or videos.



SonoMeasure

SonoMeasure features built-in measurement reference ranges for animals of different sizes, automatically alerting to potential clinical symptoms for early detection and timely intervention.



SonoMagnifier

The measurement magnifier allows you to see all the details with precision.



SonoReport

Sono-Report is a built-in veterinary report template covering ABD, Cardiac, Uro, Ophthalmic, and Repro applications, simplifying routine companion animal diagnoses for veterinarians.




Guardian of Companions Guidance for Care



Canine Kidney




Canine Biliary Sludge



TDI PW Mode

Tissue Doppler Imaging for precise assessment of myocardial motion and velocity, aiding in evaluation of heart function.



μ -Scan


Delicately engineered to reduce speckles, improve image uniformity, and enhance border continuity, providing authentic presentation of details and enhanced lesion display.




Feel Every Heartbeat Empower Every Stride



Equine Hoof Socket




Equine Tendon



Micro F

Micro F mode clearly visualizes microvascular flow around the equine embryonic vesicle, revealing fine vessel details with outstanding microflow sensitivity.



Panaromic

Real-time 2D and Color/Power Doppler Panaromic Imaging enables to extend the field of view and any mistakes can also be easily back tracked and corrected without interrupting the scan.

Ergonomic and Durable Design for
Veterinary Ultrasound Excellence



ProPet X10
Trackball (Button)



ProPet E10
Touchpad (Waterproof)



Battery Level Display



Magnetic Charging Port



Main Unit Probe Protective Cover



Silicone Panel Film



Trolley with Adjustable
Height and Tilttable Panel



Mobile Battery with High Capacity



Dust & Splash Protection

Probe Slot Dust Cover
Main Unit Dust Cover

Designed for Veterinarians

Small Probe Slot for Veterinary Use
200+ Veterinary Bodymarks
120+ Annotations Library

Enhanced Mobility

Durable Polycarbonate Trolley Case
2-in-1 Cabin Backpack & Trolley
Dedicated Probe Storage Bag